



**Vascular**<sup>®</sup>  
therapies for living

restorer

---

CoCr peripheral  
stent system

---

THE PERIPHERAL SOLUTION

# Outstanding deliverability

Based on the successful oceanus 35 balloon platform, restorer delivers 1st class performance with:

- High pushability
- Fastest balloon deflation

Enhanced trackability due to restorer's unique, flexible strut design featuring:

- Nested Peaks
- Helicoidal linkage system

The low profile design of restorer ensures increased crossability in difficult lesions.



REACH

CROSS


DELIVER



# Extraordinary radial force with less metal

Restorer demonstrates excellent vessel support


Ratio Radial Force / Wall Thickness (MPa)



Data for 6mm diameter stents. Tests performed by iVascular. Data on file at iVascular.

## Unique design

The unique design of restorer provides superior conformability to the vessel wall.



## Lower Restenosis

Utilizing an advanced manufacturing process, restorer is designed to minimize the potential for restenosis.

- Optimized metal content for maximum performance.
- Ultrasmooth electropolished surface.



## General features

- Over the wire catheter (OTW)
- Usable catheter length: 80cm & 140cm
- Guidewire compatibility: 0,035"
- Tip Profile: 0,036"
- 6F introducer compatible



CE Mark certified product

### References Catheter length: 80 cm

Stent diameter (mm)	Stent length (mm)			
	18	28	38	58
5	SPC BC35 080 050 018	SPC BC35 080 050 028	SPC BC35 080 050 038	SPC BC35 080 050 058
6	SPC BC35 080 060 018	SPC BC35 080 060 028	SPC BC35 080 060 038	SPC BC35 080 060 058
7	SPC BC35 080 070 018	SPC BC35 080 070 028	SPC BC35 080 070 038	SPC BC35 080 070 058
8	SPC BC35 080 080 018	SPC BC35 080 080 028	SPC BC35 080 080 038	SPC BC35 080 080 058
9	SPC BC35 080 090 018	SPC BC35 080 090 028	SPC BC35 080 090 038	SPC BC35 080 090 058
10	SPC BC35 080 100 018	SPC BC35 080 100 028	SPC BC35 080 100 038	SPC BC35 080 100 058

### References Catheter length: 140 cm

Stent diameter (mm)	Stent length (mm)			
	18	28	38	58
5	SPC BC35 140 050 018	SPC BC35 140 050 028	SPC BC35 140 050 038	SPC BC35 140 050 058
6	SPC BC35 140 060 018	SPC BC35 140 060 028	SPC BC35 140 060 038	SPC BC35 140 060 058
7	SPC BC35 140 070 018	SPC BC35 140 070 028	SPC BC35 140 070 038	SPC BC35 140 070 058
8	SPC BC35 140 080 018	SPC BC35 140 080 028	SPC BC35 140 080 038	SPC BC35 140 080 058
9	SPC BC35 140 090 018	SPC BC35 140 090 028	SPC BC35 140 090 038	SPC BC35 140 090 058
10	SPC BC35 140 100 018	SPC BC35 140 100 028	SPC BC35 140 100 038	SPC BC35 140 100 058

**Vascular**<sup>®</sup>  
therapies for living

Distributed by



Manufactured by:  
LVD Biotech, S.L.  
Cami de Can Ubach, 11  
08620 Sant Vicenç dels Horts  
Barcelona - Spain  
Tel: +34 936 724 711  
[www.ivascular.es](http://www.ivascular.es)  
[info@ivascular.es](mailto:info@ivascular.es)